

# MR Imaging System

OASIS

## IHE Integration Statement

 **Hitachi Medical Corporation**  
Tokyo, Japan

E1E-BM6952-5

Copyright © Hitachi Medical Corporation. 2007 - 2009, 2011. All rights reserved.

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Hitachi Medical Corporation	OASIS	V3.0	25-Mar-2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow (SWF)	Acquisition Modality	Broad Worklist Query Patient Based Worklist Query	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Consistent Presentation of Images (CPI)	Acquisition Modality	None	
	Print Composer	None	
Basic Security (SEC)	Acquisition Modality	Encryption	
	Print Composer		
	Secure Node		
Audit Trail and Node Authentication (ATNA)	Acquisition Modality	Encryption	
	Print Composer		
	Secure Node		
Consistent Time (CT)	Time Client	None	
Key Image Note (KIN)	Acquisition Modality	None	
Portable Data for Imaging (PDI)	Portable Media Creator	None	
<b>Internet address for vendor's IHE information:</b>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>			
<b>DICOM</b>			
<b>Links to general information on the IHE</b>			
In North America: <a href="http://www.ihe.net/">www.ihe.net/</a>		In Europe: <a href="http://www.ihe-europe.org/">www.ihe-europe.org/</a>	In Japan: <a href="http://www.ihe-j.org/">www.ihe-j.org/</a>