

IHE Integration Statement			
Vendor	Product Name	Version	Date
FUJIFILM Healthcare Corporation	SCENARIA View	V1.08 or later	1st Feb. 2020 Rev. 1.0
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	Broad Worklist Query	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Portable Data for Imaging (PDI)	Portable Media Creator	None	
Radiation Exposure Monitoring (REM)	Acquisition Modality	None	
Internet address for vendor's IHE information:			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM			
Links to general information on the IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org