

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Hitachi, Ltd	SCENARIO	V1.80 or later	8th Nov. 2013 Rev. 2.0
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow (SWF)	Acquisition Modality	Broad Worklist Query	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Portable Data for Imaging(PDI)	Portable Media Creator	None	
Radiation Exposure Monitoring (REM)	Acquisition Modality	None	
<b>Internet address for vendor's IHE information:</b>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>			
<b>DICOM</b>			
<b>Links to general information on the IHE</b>			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>		In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>